

## Product datasheet for **TS412263P5**

### DPS1 (PDSS1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PDSS1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212263
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<a href="#">NM_014317</a> , <a href="#">NP_055132</a>
Synonyms:	COQ1; COQ1A; COQ10D2; DPS; hDPS1; SPS; TPRT; TPT; TPT 1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_055132</a>
Locus ID:	23590
Cytogenetics:	10p12.1
Protein Pathways:	Terpenoid backbone biosynthesis



[View online »](#)