

## Product datasheet for **TS412305P5**

### **KAL1 (ANOS1) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of ANOS1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC212305
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>Target Detection Antibodies:</b>	KAL1 (ANOS1) Mouse Monoclonal Antibody [Clone ID: OT15E10] (TA808743)
<b>ACCN:</b>	<u><a href="#">NM_000216</a></u> , <u><a href="#">NP_000207</a></u>
<b>Synonyms:</b>	ADMLX; HH1; HHA; KAL; KAL1; KALIG-1; KMS; WFDC19
<b>Storage:</b>	Room Temperature
<b>Stability:</b>	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<u><a href="#">NP_000207</a></u>
<b>Locus ID:</b>	3730
<b>Cytogenetics:</b>	Xp22.31



[View online »](#)