

Product datasheet for **TS412568P5**

hnRNP F (HNRNPF) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HNRNPF, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212568
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	hnRNP F (HNRNPF) Mouse Monoclonal Antibody [Clone ID: OT15F5] (TA500804)
ACCN:	NM_001098208 , NP_001091678
Synonyms:	HNRPF; mcs94-1; OK/SW-cl.23
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001091678
Locus ID:	3185
Cytogenetics:	10q11.21



[View online »](#)