

## Product datasheet for **TS412601P5**

### KCHIP2 (KCNIP2) CytoSection

#### Product data:

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of KCNIP2, transcript variant 5, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC212601
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>Target Detection Antibodies:</b>	KCHIP2 (KCNIP2) Mouse Monoclonal Antibody [Clone ID: OT12D2] (TA807375)
<b>ACCN:</b>	<u><a href="#">NM_173194</a></u> , <u><a href="#">NP_775286</a></u>
<b>Synonyms:</b>	KCHIP2
<b>Storage:</b>	Room Temperature
<b>Stability:</b>	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<u><a href="#">NP_775286</a></u>
<b>Locus ID:</b>	30819
<b>Cytogenetics:</b>	10q24.32
<b>Protein Families:</b>	Druggable Genome, Ion Channels: Other



[View online »](#)