

Product datasheet for **TS412637P5**

DDT CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DDT, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212637
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DDT Mouse Monoclonal Antibody [Clone ID: OT13B10] (TA505591)
ACCN:	NM_001084392 , NP_001077861
Synonyms:	D-DT; DDCT; MIF-2; MIF2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001077861
Locus ID:	1652
Cytogenetics:	22q11.23



[View online »](#)