

## Product datasheet for **TS412665P5**

### Estrogen Receptor beta (ESR2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ESR2, transcript variant b, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212665
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Estrogen Receptor beta (ESR2) Mouse Monoclonal Antibody [Clone ID: OT11A5] (TA801088)
ACCN:	<a href="#">NM_001040275</a> , <a href="#">NP_001035365</a>
Synonyms:	ER-BETA; Erb; ESR-BETA; ESRB; ESTRB; NR3A2; ODG8
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001035365</a>
Locus ID:	2100
Cytogenetics:	14q23.2-q23.3
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors



[View online »](#)