

Product datasheet for **TS412883P5**

SCYL1BP1 (GORAB) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GORAB, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212883
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SCYL1BP1 (GORAB) Mouse Monoclonal Antibody [Clone ID: OTI4F11] (TA505074)
ACCN:	NM_152281 , NP_689494
Synonyms:	GO; NTKLBP1; SCYL1BP1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_689494
Locus ID:	92344
Cytogenetics:	1q24.2



[View online »](#)