

## Product datasheet for **TS412993P5**


### FOLR2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FOLR2, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212993
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	FOLR2 Mouse Monoclonal Antibody [Clone ID: OTI4G6] (TA808026)
ACCN:	<a href="#">NM_000803</a> , <a href="#">NP_000794</a>
Synonyms:	BETA-HFR; FBP; FBP/PL-1; FOLR1; FR-BETA; FR-P3; FRbeta
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_000794</a>
Locus ID:	2350
Cytogenetics:	11q13.4
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**Product images:**

Immunohistology staining on FOLR2 overexpressed cytosection TS412993P5 with mouse anti DDK clone OTI11C3 C/N [TA180144] at 1:400 dilution 4C O/N. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone [RC212993], 5 micron sections, 40x magnification