

Product datasheet for **TS413024P5**

C20orf30 (TMEM230) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TMEM230, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213024
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	C20orf30 (TMEM230) Mouse Monoclonal Antibody [Clone ID: OTI2C8] (TA504888)
ACCN:	NM_001009924 , NP_001009924
Synonyms:	C20orf30; dJ1116H23.2.1; HSPC274
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001009924
Locus ID:	29058
Cytogenetics:	20p13-p12.3
Protein Families:	Transmembrane



[View online »](#)