

Product datasheet for **TS413098P5**

TMEM16A (ANO1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ANO1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213098
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TMEM16A (ANO1) Mouse Monoclonal Antibody [Clone ID: OT11C8] (TA805184)
ACCN:	NM_018043 , NP_060513
Synonyms:	DOG1; ORAOV2; TAOS2; TMEM16A
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060513
Locus ID:	55107
Cytogenetics:	11q13.3
Protein Families:	Transmembrane



[View online »](#)