

Product datasheet for **TS413313P5**

NFkB p100 / p52 (NFKB2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NFKB2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213313
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_002502</u> , <u>NP_002493</u>
Synonyms:	CVID10; H2TF1; LYT-10; LYT10; NF-kB2; p49/p100; p52; p100
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_002493</u>
Locus ID:	4791
Cytogenetics:	10q24.32
Protein Families:	Transcription Factors
Protein Pathways:	MAPK signaling pathway, Pathways in cancer



[View online »](#)