

Product datasheet for **TS413487P5**

NFH (NEFH) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NEFH in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213487
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	NFH (NEFH) Mouse Monoclonal Antibody [Clone ID: OTI10A7] (TA801161)
ACCN:	NM_021076 , NP_066554
Synonyms:	CMT2CC; NFH
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_066554
Locus ID:	4744
Cytogenetics:	22q12.2
Protein Families:	Druggable Genome
Protein Pathways:	Amyotrophic lateral sclerosis (ALS)



[View online »](#)