

Product datasheet for TS413639P5

WNT3A CytoSection

Product data:

Product Type: CytoSections **Description:** Transient overexpression of WNT3A in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack Species: Human **Expression Host:** HEK293T **Expression cDNA Clone** TrueORF Clone RC213639 or AA Sequence: C-MYC/DDK Tag: **Detection Antibodies:** Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100) ACCN: NM 033131, NP 149122 Storage: **Room Temperature** Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. **Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. Note: This product is for research use only and is not approved for use in humans or in clinical diagnosis. **RefSeq:** NP 149122 89780 Locus ID: **Cytogenetics:** 1q42.13 **Protein Families:** Adult stem cells, Cancer stem cells, Druggable Genome, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Secreted Protein, Stem cell relevant signaling - Wnt Signaling pathway **Protein Pathways:** Basal cell carcinoma, Hedgehog signaling pathway, Melanogenesis, Pathways in cancer, Wnt signaling pathway



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn