

Product datasheet for TS413772P5

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD117 CytoSection

Product data:

Product Type: CytoSections

Description: Transient overexpression of KIT, transcript variant 2, in HEK293T cells, FFPE control for IHC,

ICC and ISH staining, 25 slides per pack

Species: Human **Expression Host:** HEK293T

Expression cDNA Clone

or AA Sequence:

TrueORF Clone RC213772

C-MYC/DDK Tag:

Detection Antibodies: DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

Target Detection

Antibodies:

CD117 Mouse Monoclonal Antibody [Clone ID: OTI9A11] (TA801205)

ACCN: NM 001093772, NP 001087241 C-Kit; CD117; MASTC; PBT; SCFR Synonyms:

Storage: Room Temperature

Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

This product is for research use only and is not approved for use in humans or in clinical Note:

diagnosis.

RefSeq: NP 001087241

Locus ID: 3815 Cytogenetics: 4q12

Protein Families: Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Protein

Kinase, Stem cell - Pluripotency, Transmembrane







Protein Pathways:

Acute myeloid leukemia, Cytokine-cytokine receptor interaction, Endocytosis, Hematopoietic cell lineage, Melanogenesis, Pathways in cancer