

## Product datasheet for **TS413781P5**

### APPBP1 (NAE1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NAE1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213781
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	APPBP1 (NAE1) Mouse Monoclonal Antibody [Clone ID: OT11E10] (TA804384)
ACCN:	<a href="#">NM_003905</a> , <a href="#">NP_003896</a>
Synonyms:	A-116A10.1; APPBP1; HPP1; ula-1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_003896</a>
Locus ID:	8883
Cytogenetics:	16q22.1
Protein Pathways:	Alzheimer's disease



[View online »](#)