

# Product datasheet for TS413809P5

## **MCM7 CytoSection**

### **Product data:**

#### **Product Type:** CytoSections **Description:** Transient overexpression of MCM7, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack Species: Human **Expression Host:** HEK293T **Expression cDNA Clone** TrueORF Clone RC213809 or AA Sequence: C-MYC/DDK Tag: **Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) ACCN: NM 005916, NP 005907 CDC47; MCM2; P1.1-MCM3; P1CDC47; P85MCM; PNAS146; PPP1R104 Synonyms: Storage: **Room Temperature** Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. **Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. Note: This product is for research use only and is not approved for use in humans or in clinical diagnosis. NP 005907 RefSeq: Locus ID: 4176 **Cytogenetics:** 7q22.1 **Protein Families: Transcription Factors Protein Pathways:** Cell cycle, DNA replication



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn