

## Product datasheet for **TS413869**

### USP54 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of USP54 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213869
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	USP54 Mouse Monoclonal Antibody [Clone ID: OTI1G6] (TA800284)
ACCN:	<a href="#">NM_152586</a> , <a href="#">NP_689799</a>
Synonyms:	bA137L10.3; bA137L10.4; C10orf29
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_689799</a>
Locus ID:	159195
Cytogenetics:	10q22.2
Protein Families:	Protease



[View online »](#)