

## Product datasheet for **TS414376P5**

### IL1RA (IL1RN) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IL1RN, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214376
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	IL1RA (IL1RN) Mouse Monoclonal Antibody [Clone ID: OTI3B1] (TA803396)
ACCN:	<a href="#">NM_173841</a> , <a href="#">NP_776213</a>
Synonyms:	DIRA; ICIL-1RA; IL-1ra; IL-1ra3; IL-1RN; IL1F3; IL1RA; IRAP; MVCD4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_776213</a>
Locus ID:	3557
Cytogenetics:	2q14.1
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)