

Product datasheet for **TS414466P5**

SLC16A2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SLC16A2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214466
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_006517</u> , <u>NP_006508</u>
Synonyms:	AHDS; DXS128; DXS128E; MCT 7; MCT7; MCT 8; MCT8; MRX22; XPCT
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_006508</u>
Locus ID:	6567
Cytogenetics:	Xq13.2
Protein Families:	Transmembrane


[View online »](#)