

Product datasheet for **TS414618P5**

ROR alpha (RORA) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RORA, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214618
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ROR alpha (RORA) Mouse Monoclonal Antibody [Clone ID: OTI2C4] (TA803201)
ACCN:	NM_002943 , NP_002934
Synonyms:	IDDECA; NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002934
Locus ID:	6095
Cytogenetics:	15q22.2
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors



[View online »](#)