

## Product datasheet for **TS414679P5**

### DIO2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DIO2, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214679
Tag:	Tag Free
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_013989</a> , <a href="#">NP_054644</a>
Synonyms:	5DII; D2; DIOII; SELENOY; SeLY; TXDI2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_054644</a>
Locus ID:	1734
Cytogenetics:	14q31.1
Protein Families:	Druggable Genome



[View online »](#)