

Product datasheet for **TS414978P5**

KIF6 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KIF6 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214978
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	KIF6 Mouse Monoclonal Antibody [Clone ID: OTI5G2] (TA811816)
ACCN:	<u>NM_145027</u> , <u>NP_659464</u>
Synonyms:	C6orf102; dj137F1.4; dj188D3.1; dj1043E3.1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_659464</u>
Locus ID:	221458
Cytogenetics:	6p21.2
Protein Families:	Druggable Genome


[View online »](#)