

Product datasheet for **TS414999P5**

ILF3 CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of ILF3, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC214999 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | NM_012218 , NP_036350 |
| Synonyms: | CBTF; DRBF; DRBP76; MMP4; MPHOSPH4; MPP4; MPP41 10; NF-AT-90; NF90; NF90a; NF90b; NF90c; NF90ctv; NF110; NF110b; NFAR; NFAR-1; NFAR-2; NFAR2; NFAR90; NFAR110; TCP80; TCP110 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_036350 |
| Locus ID: | 3609 |
| Cytogenetics: | 19p13.2 |
| Protein Families: | Druggable Genome, Transcription Factors |



[View online »](#)