

Product datasheet for **TS415028P5**

LYK5 (STRADA) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of STRADA, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215028
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001003788 , NP_001003788
Synonyms:	LYK5; NY-BR-96; PMSE; Stlk; STRAD; STRAD alpha
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001003788
Locus ID:	92335
Cytogenetics:	17q23.3
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	mTOR signaling pathway



[View online »](#)