

## Product datasheet for **TS415105P5**

### GSG1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GSG1, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215105
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GSG1 Mouse Monoclonal Antibody [Clone ID: OT11H9] (TA812001)
ACCN:	<u><a href="#">NM_153823</a></u> , <u><a href="#">NP_722545</a></u>
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_722545</a></u>
Locus ID:	83445
Cytogenetics:	12p13.1
Protein Families:	Transmembrane



[View online »](#)