

Product datasheet for **TS415126P5**

Transglutaminase 5 (TGM5) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TGM5, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215126
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_201631 , NP_963925
Synonyms:	PSS2; TGASE5; TGASEX; TGM6; TGMX; TGX
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_963925
Locus ID:	9333
Cytogenetics:	15q15.2
Protein Families:	Druggable Genome



[View online »](#)