

## Product datasheet for **TS415149P5**

### H4C4 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of H4C4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215149
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_003539</a></u> , <u><a href="#">NP_003530</a></u>
Synonyms:	dj221C16.9; H4-16; H4/b; H4C1; H4C2; H4C3; H4C5; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C14; H4C15; H4FB; HIST1H4D
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_003530</a></u>
Locus ID:	8360
Cytogenetics:	6p22.2
Protein Pathways:	Systemic lupus erythematosus



[View online »](#)