

## Product datasheet for **TS415204P5**

### DPP8 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DPP8, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215204
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DPP8 Mouse Monoclonal Antibody [Clone ID: OTI1D2] (TA504366)
ACCN:	<a href="#">NM_017743</a> , <a href="#">NP_060213</a>
Synonyms:	DP8; DPRP-1; DPRP1; MST097; MSTP097; MSTP135; MSTP141
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_060213</a>
Locus ID:	54878
Cytogenetics:	15q22.31
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)