

Product datasheet for **TS415286P5**

Lipin 3 (LPIN3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LPIN3 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215286
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Lipin 3 (LPIN3) Mouse Monoclonal Antibody [Clone ID: OTI1D7] (TA804944)
ACCN:	<u>NM_022896</u> , <u>NP_075047</u>
Synonyms:	dj620E11.2; LIPN3L; SMP2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_075047</u>
Locus ID:	64900
Cytogenetics:	20q12



[View online »](#)