

## Product datasheet for **TS415291**

### **B7H3 (CD276) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of CD276, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC215291
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>Target Detection Antibodies:</b>	B7H3 (CD276) Mouse Monoclonal Antibody [Clone ID: OT11A10] (TA813302)
<b>ACCN:</b>	<a href="#">NM_025240</a> , <a href="#">NP_079516</a>
<b>Synonyms:</b>	4Ig-B7-H3; B7-H3; B7H3; B7RP-2
<b>Storage:</b>	Room Temperature
<b>Stability:</b>	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<a href="#">NP_079516</a>
<b>Locus ID:</b>	80381
<b>Cytogenetics:</b>	15q24.1
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Cell adhesion molecules (CAMs)



[View online »](#)