

Product datasheet for **TS415300**

Ceramide synthase 2 (CERS2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CERS2, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215300
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Ceramide synthase 2 (CERS2) Mouse Monoclonal Antibody [Clone ID: OTI3D9] (TA809918)
ACCN:	NM_181746 , NP_859530
Synonyms:	L3; LASS2; SP260; TMSG1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_859530
Locus ID:	29956
Cytogenetics:	1q21.3
Protein Families:	Transcription Factors, Transmembrane



[View online »](#)