

Product datasheet for **TS415356P5**

KRTAP2-4 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KRTAP2-4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215356
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	KRTAP2-4 Mouse Monoclonal Antibody [Clone ID: OT13B4] (TA503766)
ACCN:	<u>NM_033184, NP_149440</u>
Synonyms:	KAP2.1B; KAP2.3; KAP2.4; KRTAP2-3; KRTAP2.3; KRTAP2.4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_149440</u>
Locus ID:	85294
Cytogenetics:	17q21.2


[View online »](#)