

Product datasheet for **TS415488P5**

mtTFA (TFAM) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TFAM in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215488
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_003201 , NP_003192
Synonyms:	MTDPS15; MTTF1; MTTFA; TCF6; TCF6L1; TCF6L2; TCF6L3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003192
Locus ID:	7019
Cytogenetics:	10q21.1
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Huntington's disease



[View online »](#)