

Product datasheet for **TS415596P5**

H4C12 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of H4C12 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215596
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_003541</u> , <u>NP_003532</u>
Synonyms:	dj160A22.1; H4-16; H4/d; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4C11; H4C13; H4C14; H4C15; H4F2iii; H4FD; HIST1H4K
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_003532</u>
Locus ID:	8362
Cytogenetics:	6p22.1
Protein Pathways:	Systemic lupus erythematosus



[View online »](#)