

Product datasheet for TS415661P5

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MDA5 (IFIH1) CytoSection

Product data:

Product Type: CytoSections

Description: Transient overexpression of IFIH1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining,

25 slides per pack

Species: Human
Expression Host: HEK293T

Expression cDNA Clone

or AA Sequence:

TrueORF Clone RC215661

Tag: C-MYC/DDK

Detection Antibodies: Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)

Target Detection

Antibodies:

MDA5 (IFIH1) Mouse Monoclonal Antibody [Clone ID: OTI8C10] (TA803667)

ACCN: <u>NM 022168</u>, <u>NP 071451</u>

Synonyms: AGS7; Hlcd; IDDM19; MDA-5; MDA5; RLR-2; SGMRT1

Storage: Room Temperature

Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

Note: This product is for research use only and is not approved for use in humans or in clinical

diagnosis.

RefSeq: NP 071451

Locus ID: 64135 Cytogenetics: 2q24.2

Protein Pathways: RIG-I-like receptor signaling pathway

