

Product datasheet for **TS415726P5**

HMG2L1 (HMGXB4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HMGXB4, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215726
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HMG2L1 (HMGXB4) Mouse Monoclonal Antibody [Clone ID: OTI1H9] (TA801840)
ACCN:	NM_014250 , NP_055065
Synonyms:	high-mobility group protein 2-like 1; HMG2L1; HMGBCG; HMG box domain containing 4; OTTHUMP00000028778; THC211630
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_055065
Locus ID:	10042
Cytogenetics:	22q12.3
Protein Families:	Transcription Factors



[View online »](#)