

## Product datasheet for **TS415953**

### HIF3 alpha (HIF3A) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HIF3A, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215953
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HIF3 alpha (HIF3A) Mouse Monoclonal Antibody [Clone ID: OTI2D2] (TA800720)
ACCN:	<a href="#">NM_022462</a> , <a href="#">NP_071907</a>
Synonyms:	bHLHe17; HIF-3A; HIF3-alpha-1; IPAS; MOP7; PASD7
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_071907</a>
Locus ID:	64344
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)