

Product datasheet for **TS416051P5**

ZNF197 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZNF197, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216051
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001024855 , NP_001020026
Synonyms:	D3S1363E; P18; VHLaK; ZKSCAN9; ZNF20; ZNF166; ZSCAN41
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001020026
Locus ID:	10168
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)