

## Product datasheet for **TS416146P5**

### KCNG4 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCNG4, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216146
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_133490</a></u> , <u><a href="#">NP_597997</a></u>
Synonyms:	KV6.3; KV6.3, KV6.4, MGC4558, MGC129609; KV6.4; MGC4558; MGC129609; potassium voltage-gated channel, subfamily G, member 4; potassium voltage-gated channel, subfamily G, member 4, isoform 1; voltage-gated potassium channel Kv6.3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_597997</a></u>
Locus ID:	93107
Cytogenetics:	16q24.1
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane



[View online »](#)