

## Product datasheet for **TS416176P5**

### PPAR alpha (PPARA) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PPARA, transcript variant 5, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216176
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PPAR alpha (PPARA) Mouse Monoclonal Antibody [Clone ID: OT11E8] (TA503540)
ACCN:	<a href="#">NM_005036</a> , <a href="#">NP_005027</a>
Synonyms:	hPPAR; NR1C1; PPAR; PPAR-alpha; PPARalpha
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_005027</a>
Locus ID:	5465
Cytogenetics:	22q13.31
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Protein Pathways:	Adipocytokine signaling pathway, PPAR signaling pathway



[View online »](#)