

## Product datasheet for **TS416204P5**

### TBX1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TBX1, transcript variant C, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216204
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TBX1 Mouse Monoclonal Antibody [Clone ID: OTI1C2] (TA507326)
ACCN:	<a href="#">NM_080647</a> , <a href="#">NP_542378</a>
Synonyms:	CAFS; CATCH22; CTHM; DGCR; DGS; DORV; TBX1C; TGA; VCF; VCFS
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_542378</a>
Locus ID:	6899
Cytogenetics:	22q11.21
Protein Families:	Transcription Factors



[View online »](#)