

Product datasheet for **TS416263P5**

TXNRD1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TXNRD1, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216263
Tag:	Tag Free
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TXNRD1 Mouse Monoclonal Antibody [Clone ID: OT13E1] (TA811364)
ACCN:	NM_182743 , NP_877420
Synonyms:	GRIM-12; TR; TR1; TRXR1; TXNR
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_877420
Locus ID:	7296
Cytogenetics:	12q23.3
Protein Families:	Druggable Genome
Protein Pathways:	Pyrimidine metabolism



[View online »](#)