

Product datasheet for **TS416801P5**

UBE1C (UBA3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UBA3, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216801
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	UBE1C (UBA3) Mouse Monoclonal Antibody [Clone ID: OTI7H6] (TA811903)
ACCN:	NM_198195 , NP_937838
Synonyms:	hUBA3; NAE2; UBE1C
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_937838
Locus ID:	9039
Cytogenetics:	3p14.1
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)