

## Product datasheet for **TS416836P5**

### PAX7 CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of PAX7, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack  |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC216836  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| Target Detection Antibodies:          | PAX7 Mouse Monoclonal Antibody [Clone ID: OTI4A7] (TA811661)  |
| ACCN:                                 | <a href="#">NM_013945</a> , <a href="#">NP_039236</a>   |
| Synonyms:                             | HUP1; MYOSCO; PAX7B; RMS2   |
| Storage:                              | Room Temperature  |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <a href="#">NP_039236</a>   |
| Locus ID:                             | 5081  |
| Cytogenetics:                         | 1p36.13   |
| Protein Families:                     | Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Transcription Factors  |



[View online »](#)