

## Product datasheet for **TS416869P5**

### NUMBL CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NUMBL in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216869
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	NUMBL Mouse Monoclonal Antibody [Clone ID: OTI4A1] (TA812189)
ACCN:	<a href="#">NM_004756</a> , <a href="#">NP_004747</a>
Synonyms:	CAG3A; CTG3a; NBL; NUMB-R; NUMBLIKE; NUMBR; TNRC23
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_004747</a>
Locus ID:	9253
Cytogenetics:	19q13.2
Protein Pathways:	Notch signaling pathway



[View online »](#)