

## Product datasheet for **TS416890P5**

### **C14orf166B (LRRC74A) CytoSection**

#### **Product data:**

|  |   |
|--|---|
| <b>Product Type:</b>                         | CytoSections  |
| <b>Description:</b>                          | Transient overexpression of C14orf166B in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack   |
| <b>Species:</b>                              | Human   |
| <b>Expression Host:</b>                      | HEK293T   |
| <b>Expression cDNA Clone or AA Sequence:</b> | TrueORF Clone RC216890  |
| <b>Tag:</b>                                  | C-MYC/DDK   |
| <b>Detection Antibodies:</b>                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| <b>ACCN:</b>                                 | <u><a href="#">NM_001081002</a></u>   |
| <b>Synonyms:</b>                             | chromosome 14 open reading frame 166B   |
| <b>Storage:</b>                              | Room Temperature  |
| <b>Stability:</b>                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| <b>Preparation:</b>                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| <b>Note:</b>                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| <b>Locus ID:</b>                             | 145497  |
| <b>Cytogenetics:</b>                         | 14q24.3   |



[View online »](#)