

Product datasheet for **TS416929P5**

SCD5 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SCD5, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216929
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_024906 , NP_079182
Synonyms:	ACOD4; FADS4; HSCD5; SCD2; SCD4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_079182
Locus ID:	79966
Cytogenetics:	4q21.22
Protein Families:	Transmembrane
Protein Pathways:	Biosynthesis of unsaturated fatty acids, PPAR signaling pathway



[View online »](#)