

## Product datasheet for **TS416947P5**


### IL13 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IL13 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216947
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	IL13 Mouse Monoclonal Antibody [Clone ID: OTI6D3] (TA804004)
ACCN:	<u>NM_002188, NP_002179</u>
Synonyms:	IL-13; P600
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_002179</u>
Locus ID:	3596
Cytogenetics:	5q31.1
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	Asthma, Cytokine-cytokine receptor interaction, Fc epsilon RI signaling pathway, Jak-STAT signaling pathway



[View online »](#)

**Product images:**

Immunohistology staining on IL13 overexpressed cytosection TS416947P5 with rabbit anti DDK clone OTIR5G2 C/N [TA592569] at 1:2000 dilution 5m RT. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone , 5 micron sections, 40x magnification