

## Product datasheet for **TS417121**

### FGF13 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FGF13, transcript variant 6, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217121
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	FGF13 Mouse Monoclonal Antibody [Clone ID: OTI7E8] (TA807761)
ACCN:	<a href="#">NM_033642</a> , <a href="#">NP_378668</a>
Synonyms:	DEE90; FGF-13; FGF2; FHF-2; FHF2; LINC00889
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_378668</a>
Locus ID:	2258
Cytogenetics:	Xq26.3-q27.1
Protein Families:	Secreted Protein
Protein Pathways:	MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton



[View online »](#)