

Product datasheet for **TS417297P5**

CPT1B CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of CPT1B, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC217297 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | CPT1B Mouse Monoclonal Antibody [Clone ID: OTI2A6] (TA810597) |
| ACCN: | NM_004377 , NP_004368 |
| Synonyms: | CPT1-M; CPT1M; CPTI; CPTI-M; M-CPT1; MCCPT1; MCPT1 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_004368 |
| Locus ID: | 1375 |
| Cytogenetics: | 22q13.33 |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | Adipocytokine signaling pathway, Fatty acid metabolism, PPAR signaling pathway |



[View online »](#)